

IBP4 Polyclonal Antibody

Cat No: HR1AP10265

For research use only

Overview

Product Name	IBP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	IGFBP4 IBP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>insulin like growth factor binding protein 4(IGFBP4) Homo sapiens This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma in both glycosylated and non-glycosylated forms. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors. [provided by RefSeq, Jul 2008],</p>