

## GTSE1 Polyclonal Antibody

Cat No: HR1AP10261

For research use only

### Overview

Product Name	GTSE1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	GTSE1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G2 and S-phase expressed 1(GTSE1) Homo sapiens The protein encoded by this gene is only expressed in the S and G2 phases of the cell cycle, where it colocalizes with cytoplasmic tubulin and microtubules. In response to DNA damage, the encoded protein accumulates in the nucleus and binds the tumor suppressor protein p53, shuttling it out of the nucleus and repressing its ability to induce apoptosis. [provided by RefSeq, Jul 2008],