

# FGFP3 Polyclonal Antibody

Cat No: HR1AP10257

For research use only

## Overview

Product Name	FGFP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	FGFBP3 C10orf13 PSEC0101
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the fibroblast growth factor-binding protein family.,subunit:Interacts with FGF1.,