

FRAT1 Polyclonal Antibody

Cat No: HR1AP10259

For research use only

Overview

Product Name	FRAT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	FRAT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	frequently rearranged in advanced T-cell lymphomas 1(FRAT1) Homo sapiens The protein encoded by this gene belongs to the GSK-3-binding protein family. The protein inhibits GSK-3-mediated phosphorylation of beta-catenin and positively regulates the Wnt signaling pathway. It may function in tumor progression and in lymphomagenesis. [provided by RefSeq, Oct 2008],