

FGF21 Polyclonal Antibody

Cat No: HR1AP10256

For research use only

Overview

Product Name	FGF21 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	FGF21 UNQ3115/PRO10196
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>fibroblast growth factor 21(FGF21) Homo sapiens This gene encodes a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities and are involved in a variety of biological processes. This protein is a secreted endocrine factor that functions as a major metabolic regulator. The encoded protein stimulates the uptake of glucose in adipose tissue. [provided by RefSeq, Mar 2016],</p>