

EP15R Polyclonal Antibody

Cat No: HR1AP10252

For research use only

Overview

Product Name	EP15R Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	EPS15L1 EPS15R
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Seems to be a constitutive component of clathrin-coated pits that is required for receptor-mediated endocytosis.,PTM:Phosphorylated on tyrosine residues by EGFR.,similarity:Contains 1 EF-hand domain.,similarity:Contains 3 EH domains.,subcellular location:Localized to plasma membrane coated pits.,subunit:Interacts with EPS15, AGFG1/HRB and AGFG2/HRBL. Associates with the clathrin-associated adapter protein complex 2 (AP-2).</p>