

CALY Polyclonal Antibody

Cat No: HR1AP10250

For research use only

Overview

Product Name	CALY Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	CALY DRD1IP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	calcyon neuron specific vesicular protein(CALY) Homo sapiens The protein encoded by this gene is a type II single transmembrane protein. It is required for maximal stimulated calcium release after stimulation of purinergic or muscarinic but not beta-adrenergic receptors. The encoded protein interacts with D1 dopamine receptor and may interact with other DA receptor subtypes and/or GPCRs. [provided by RefSeq, Jul 2008],