

## DNJB5 Polyclonal Antibody

Cat No: HR1AP10251

For research use only

### Overview

Product Name	DNJB5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	DNAJB5 HSC40
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DnaJ heat shock protein family (Hsp40) member B5(DNAJB5) Homo sapiens DNAJB5 belongs to the evolutionarily conserved DNAJ/HSP40 protein family. For background information on the DNAJ family, see MIM 608375.[supplied by OMIM, Mar 2008],