

## CRLF1 Polyclonal Antibody

Cat No: HR1AP10244

For research use only

### Overview

Product Name	CRLF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	CRLF1 UNQ288/PRO327
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cytokine receptor like factor 1(CRLF1) Homo sapiens This gene encodes a member of the cytokine type I receptor family. The protein forms a secreted complex with cardiotrophin-like cytokine factor 1 and acts on cells expressing ciliary neurotrophic factor receptors. The complex can promote survival of neuronal cells. Mutations in this gene result in Crisponi syndrome and cold-induced sweating syndrome. [provided by RefSeq, Oct 2009],</p>