

## CD320 Polyclonal Antibody

Cat No: HR1AP10228

For research use only

### Overview

Product Name	CD320 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	CD320 8D6A UNQ198/PRO224
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>CD320 molecule(CD320) Homo sapiens This gene encodes the transcobalamin receptor that is expressed at the cell surface. It mediates the cellular uptake of transcobalamin bound cobalamin (vitamin B12), and is involved in B-cell proliferation and immunoglobulin secretion. Mutations in this gene are associated with methylmalonic aciduria. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2011],</p>