

## CO8A1 Polyclonal Antibody

Cat No: HR1AP10224

For research use only

### Overview

Product Name	CO8A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	COL8A1 C3orf7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	collagen type VIII alpha 1 chain(COL8A1) Homo sapiens This gene encodes one of the two alpha chains of type VIII collagen. The gene product is a short chain collagen and a major component of the basement membrane of the corneal endothelium. The type VIII collagen fibril can be either a homo- or a heterotrimer. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by RefSeq, Dec 2011],