

MAFB Polyclonal Antibody

Cat No: HR1AP10221

For research use only

Overview

Product Name	MAFB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	MAFB KRML
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MAF bZIP transcription factor B(MAFB) Homo sapiens The protein encoded by this gene is a basic leucine zipper (bZIP) transcription factor that plays an important role in the regulation of lineage-specific hematopoiesis. The encoded nuclear protein represses ETS1-mediated transcription of erythroid-specific genes in myeloid cells. This gene contains no introns. [provided by RefSeq, Jul 2008],