

WNK2 Polyclonal Antibody

Cat No: HR1AP10220

For research use only

Overview

Product Name	WNK2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	252kDa
GeneID?Human?	WNK2 KIAA1760 PRKWNK2 SDCCAG43 P/OKcl.13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	WNK lysine deficient protein kinase 2(WNK2) Homo sapiens The protein encoded by this gene is a cytoplasmic serine-threonine kinase that belongs to the protein kinase superfamily. The protein plays an important role in the regulation of electrolyte homeostasis, cell signaling survival, and proliferation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013],