

# TR150 Polyclonal Antibody

Cat No: HR1AP10219

For research use only

## Overview

Product Name	TR150 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	105kDa
GeneID?Human?	THRAP3 TRAP150
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Plays a role in transcriptional coactivation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Subunit of the large multiprotein complex TRAP.,tissue specificity:Ubiquitous.,