

## TRI35 Polyclonal Antibody

Cat No: HR1AP10215

For research use only

### Overview

Product Name	TRI35 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	TRIM35 HLS5 KIAA1098
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tripartite motif containing 35(TRIM35) Homo sapiens The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The function of this protein has not been identified. [provided by RefSeq, Jul 2008],