

## RBPJL Polyclonal Antibody

Cat No: HR1AP10211

For research use only

### Overview

Product Name	RBPJL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	RBPJL RBPL RBPSUHL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>recombination signal binding protein for immunoglobulin kappa J region like(RBPJL) Homo sapiens This gene encodes a member of the suppressor of hairless protein family. A similar protein in mouse is a transcription factor that binds to DNA sequences almost identical to that bound by the Notch receptor signaling pathway transcription factor recombining binding protein J. The mouse protein has been shown to activate transcription in concert with Epstein-Barr virus nuclear antigen-2. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],</p>