

MBIP1 Polyclonal Antibody

Cat No: HR1AP10207

For research use only

Overview

Product Name	MBIP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	MBIP BM-015
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Inhibits the MAP3K12 activity to induce the activation of the JNK/SAPK pathway.,subcellular location:Shows a cytoplasmic localization when co-expressed with MAP3K12.,tissue specificity:Ubiquitous. High expression seen in the heart and lung.,