

BARH2 Polyclonal Antibody

Cat No: HR1AP10209

For research use only

Overview

Product Name	BARH2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	BARHL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Potential regulator of neural basic helix-loop-helix genes.,similarity:Belongs to the BAR homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,