

ABC3C Polyclonal Antibody

Cat No: HR1AP10206

For research use only

Overview

Product Name	ABC3C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	APOBEC3C APOBEC1L PBI
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>apolipoprotein B mRNA editing enzyme catalytic subunit 3C(APOBEC3C) Homo sapiens This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control. [provided by RefSeq, Jul 2008].</p>