

ZMAT3 Polyclonal Antibody

Cat No: HR1AP10204

For research use only

Overview

Product Name	ZMAT3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	ZMAT3 PAG608 WIG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>zinc finger matrin-type 3(ZMAT3) Homo sapiens This gene encodes a protein containing three zinc finger domains and a nuclear localization signal. The mRNA and the protein of this gene are upregulated by wildtype p53 and overexpression of this gene inhibits tumor cell growth, suggesting that this gene may have a role in the p53-dependent growth regulatory pathway. Alternative splicing of this gene results in two transcript variants encoding two isoforms differing in only one amino acid. [provided by RefSeq, Jul 2008],</p>