

## KPSH2 Polyclonal Antibody

Cat No: HR1AP10199

For research use only

### Overview

Product Name	KPSH2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	PSKH2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,caution:Asn-183 is present instead of the conserved Asp which is expected to be an active site residue.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,