

# TF2H4 Polyclonal Antibody

Cat No: HR1AP10196

For research use only

## Overview

Product Name	TF2H4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	GTF2H4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Component of the core-TFIIF basal transcription factor involved in nucleotide excision repair (NER) of DNA and, when complexed to CAK, in RNA transcription by RNA polymerase II.,similarity:Belongs to the TFB2 family.,subunit:One of the six subunits forming the core-TFIIF basal transcription factor.,