

LMO7 Polyclonal Antibody

Cat No: HR1AP10192

For research use only

Overview

Product Name	LMO7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	185kDa
GeneID?Human?	LMO7 FBX20 FBXO20 KIAA0858
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	LIM domain 7(LMO7) Homo sapiens This gene encodes a protein containing a calponin homology (CH) domain, a PDZ domain, and a LIM domain, and may be involved in protein-protein interactions. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene, however, the full-length nature of some variants is not known. [provided by RefSeq, Jan 2009],