

FAKD5 Polyclonal Antibody

Cat No: HR1AP10186

For research use only

Overview

Product Name	FAKD5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	FASTKD5 KIAA1792
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,