

# M3K15 Polyclonal Antibody

Cat No: HR1AP10184

For research use only

## Overview

Product Name	M3K15 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	144kDa
GeneID?Human?	MAP3K15 ASK3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>mitogen-activated protein kinase kinase kinase 15(MAP3K15) Homo sapiens The protein encoded by this gene is a member of the mitogen-activated protein kinase (MAPK) family. These family members function in a protein kinase signal transduction cascade, where an activated MAPK kinase kinase (MAP3K) phosphorylates and activates a specific MAPK kinase (MAP2K), which then activates a specific MAPK. This MAP3K protein plays an essential role in apoptotic cell death triggered by cellular stresses. [provided by RefSeq, Jul 2010],</p>