

## PXDC2 Polyclonal Antibody

Cat No: HR1AP10182

For research use only

### Overview

Product Name	PXDC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GeneID?Human?	PLXDC2 TEM7R UNQ2514/PRO6003
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in tumor angiogenesis.,similarity:Belongs to the plexin family.,similarity:Contains 1 PSI domain.,subunit:Interacts with CTTN.,tissue specificity:Expressed in the endothelial cells of the stroma but not in the endothelial cells of normal colonic tissue.,