



KKLC1 Polyclonal Antibody

Cat No: HR1AP10180

For research use only

Overview

Product Name	KKLC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	KKLC1 CXorf61
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tissue specificity:Specifically expressed in testis. Expressed by cancer cell lines.,