

# TAL2 Polyclonal Antibody

Cat No: HR1AP10177

For research use only

## Overview

Product Name	TAL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	TAL2 BHLHA19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TAL bHLH transcription factor 2(TAL2) Homo sapiens This intronless gene encodes a helix-loop-helix protein. Translocations between this gene on chromosome 9 and the T-cell receptor beta-chain locus on chromosome 7 have been associated with activation of the T-cell acute lymphocytic leukemia 2 gene and T-cell acute lymphoblastic leukemia. [provided by RefSeq, Mar 2009],