

# PGFRL Polyclonal Antibody

Cat No: HR1AP10174

For research use only

## Overview

Product Name	PGFRL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	PDGFRL PRLTS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	platelet derived growth factor receptor like(PDGFRL) Homo sapiens This gene encodes a protein with significant sequence similarity to the ligand binding domain of platelet-derived growth factor receptor beta. Mutations in this gene, or deletion of a chromosomal segment containing this gene, are associated with sporadic hepatocellular carcinomas, colorectal cancers, and non-small cell lung cancers. This suggests this gene product may function as a tumor suppressor. [provided by RefSeq, Jul 2008],