

# NOCT Polyclonal Antibody

Cat No: HR1AP10172

For research use only

## Overview

Product Name	NOCT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	CCRN4L CCR4 NOC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>nocturnin(NOCT) Homo sapiens The protein encoded by this gene is highly similar to Nocturnin, a gene identified as a circadian clock regulated gene in <i>Xenopus laevis</i>. This protein and Nocturnin protein share similarity with the C-terminal domain of a yeast transcription factor, carbon catabolite repression 4 (CCR4). The mRNA abundance of a similar gene in mouse has been shown to exhibit circadian rhythmicity, which suggests a role for this protein in clock function or as a circadian clock effector. [provided by RefSeq, Jul 2008],</p>