

## DNJC3 Polyclonal Antibody

Cat No: HR1AP10168

For research use only

### Overview

Product Name	DNJC3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	DNAJC3 P58IPK PRKRI
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DnaJ heat shock protein family (Hsp40) member C3(DNAJC3) Homo sapiens This gene encodes a protein with multiple tetratricopeptide repeat (TPR) motifs as well as the highly conserved J domain found in DNAJ chaperone family members. It is a member of the tetratricopeptide repeat family of proteins and acts as an inhibitor of the interferon-induced, dsRNA-activated protein kinase (PKR). [provided by RefSeq, Jul 2010],