

## ZBT17 Polyclonal Antibody

Cat No: HR1AP10167

For research use only

### Overview

Product Name	ZBT17 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	88kDa
GeneID?Human?	ZBTB17 MIZ1 ZNF151 ZNF60
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger and BTB domain containing 17(ZBTB17) Homo sapiens This gene encodes a zinc finger protein involved in the regulation of c-myc. The symbol MIZ1 has also been associated with PIAS2 which is a different gene located on chromosome 18. [provided by RefSeq, Jul 2008],