

WIT1 Polyclonal Antibody

Cat No: HR1AP10160

For research use only

Overview

Product Name	WIT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	WT1-AS WIT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>WT1 antisense RNA(WT1-AS) Homo sapiens This gene is located upstream of the Wilms tumor 1 (WT1) gene; these two genes are bi-directionally transcribed from the same promoter region. This gene is imprinted in kidney, with preferential expression from the paternal allele. Imprinting defects at chromosome 11p13 may contribute to tumorigenesis. [provided by RefSeq, May 2014],</p>