

# FOXM1 Polyclonal Antibody

Cat No: HR1AP10161

For research use only

## Overview

Product Name	FOXM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	FOXM1 FKHL16 HFH11 MPP2 WIN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	forkhead box M1(FOXM1) Homo sapiens The protein encoded by this gene is a transcriptional activator involved in cell proliferation. The encoded protein is phosphorylated in M phase and regulates the expression of several cell cycle genes, such as cyclin B1 and cyclin D1. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011],