

## BNC1 Polyclonal Antibody

Cat No: HR1AP10159

For research use only

### Overview

Product Name	BNC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	109kDa
GeneID?Human?	BNC1 BNC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>basonuclin 1(BNC1) Homo sapiens The protein encoded by this gene is a zinc finger protein present in the basal cell layer of the epidermis and in hair follicles. It is also found in abundance in the germ cells of testis and ovary. This protein is thought to play a regulatory role in keratinocyte proliferation and it may also be a regulator for rRNA transcription. Alternative splicing of this gene results in multiple transcript variants, and multiple polyadenylation sites are indicated.[provided by RefSeq, Jul 2014],</p>