

## U2AF1 Polyclonal Antibody

Cat No: HR1AP10158

For research use only

### Overview

Product Name	U2AF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	U2AF1 U2AF35 U2AFBP FP793
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>U2 small nuclear RNA auxiliary factor 1 (U2AF1) Homo sapiens This gene belongs to the splicing factor SR family of genes. U2 auxiliary factor, comprising a large and a small subunit, is a non-snRNP protein required for the binding of U2 snRNP to the pre-mRNA branch site. This gene encodes the small subunit which plays a critical role in both constitutive and enhancer-dependent RNA splicing by directly mediating interactions between the large subunit and proteins bound to the enhancers. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],</p>