

## SPT4H Polyclonal Antibody

Cat No: HR1AP10156

For research use only

### Overview

Product Name	SPT4H Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	SUPT4H1 SPT4H SUPT4H
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SPT4 homolog, DSIF elongation factor subunit(SUPT4H1) Homo sapiens This gene encodes the small subunit of DRB (5,6-dichloro-1-beta-d-ribofuranosylbenzimidazole) sensitivity-inducing factor (DSIF) complex, which regulates mRNA processing and transcription elongation by RNA polymerase II. The encoded protein is localized to the nucleus and interacts with the large subunit (SUPT5H) to form the DSIF complex. Related pseudogenes have been identified on chromosomes 2 and 12. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Nov 2012],</p>