

SMD3 Polyclonal Antibody

Cat No: HR1AP10155

For research use only

Overview

Product Name	SMD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	SNRPD3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	small nuclear ribonucleoprotein D3 polypeptide(SNRPD3) Homo sapiens This gene encodes a core component of the spliceosome, which is a nuclear ribonucleoprotein complex that functions in pre-mRNA splicing. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],