

## VSX2 Polyclonal Antibody

Cat No: HR1AP10154

For research use only

### Overview

Product Name	VSX2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	VSX2 CHX10 HOX10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	visual system homeobox 2(VSX2) Homo sapiens This gene encodes a homeobox protein originally described as a retina-specific transcription factor. Mutations in this gene are associated with microphthalmia, cataracts and iris abnormalities. [provided by RefSeq, Oct 2009],