

## ID1 Polyclonal Antibody

Cat No: HR1AP10150

For research use only

### Overview

Product Name	ID1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	ID1 BHLHB24 ID
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>inhibitor of DNA binding 1, HLH protein(ID1) Homo sapiens The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with members of the basic HLH family of transcription factors. The encoded protein has no DNA binding activity and therefore can inhibit the DNA binding and transcriptional activation ability of basic HLH proteins with which it interacts. This protein may play a role in cell growth, senescence, and differentiation. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>