

# MK04 Polyclonal Antibody

Cat No: HR1AP10148

For research use only

## Overview

Product Name	MK04 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	MAPK4 ERK4 PRKM4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mitogen-activated protein kinase 4(MAPK4) Homo sapiens Mitogen-activated protein kinase 4 is a member of the mitogen-activated protein kinase family. Tyrosine kinase growth factor receptors activate mitogen-activated protein kinases which then translocate into the nucleus and phosphorylate nuclear targets. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014],