

## EGR2 Polyclonal Antibody

Cat No: HR1AP10145

For research use only

### Overview

Product Name	EGR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	EGR2 KROX20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	early growth response 2(EGR2) Homo sapiens The protein encoded by this gene is a transcription factor with three tandem C2H2-type zinc fingers. Defects in this gene are associated with Charcot-Marie-Tooth disease type 1D (CMT1D), Charcot-Marie-Tooth disease type 4E (CMT4E), and with Dejerine-Sottas syndrome (DSS). Multiple transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Oct 2008],