

CDKL5 Polyclonal Antibody

Cat No: HR1AP10141

For research use only

Overview

Product Name	CDKL5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	113kDa
GeneID?Human?	CDKL5 STK9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cyclin dependent kinase like 5(CDKL5) Homo sapiens This gene is a member of Ser/Thr protein kinase family and encodes a phosphorylated protein with protein kinase activity. Mutations in this gene have been associated with X-linked infantile spasm syndrome (ISSX), also known as X-linked West syndrome, and Rett syndrome (RTT). Alternate transcriptional splice variants have been characterized. [provided by RefSeq, Jul 2008],</p>