

PITX3 Polyclonal Antibody

Cat No: HR1AP10139

For research use only

Overview

Product Name	PITX3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	PITX3 PTX3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>paired like homeodomain 3(PITX3) Homo sapiens This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. Members of this family act as transcription factors. This protein is involved in lens formation during eye development. Mutations of this gene have been associated with anterior segment mesenchymal dysgenesis and congenital cataracts. [provided by RefSeq, Jul 2008],</p>