

# NKRF Polyclonal Antibody

Cat No: HR1AP10134

For research use only

## Overview

Product Name	NKRF Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	NKRF ITBA4 NRF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NFKB repressing factor(NKRF) Homo sapiens This gene encodes a transcriptional repressor that interacts with specific negative regulatory elements to mediate transcriptional repression of certain nuclear factor kappa B responsive genes. The protein localizes predominantly to the nucleolus with a small fraction found in the nucleoplasm and cytoplasm. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010],