



# LSM1 Polyclonal Antibody

Cat No: HR1AP10133

For research use only

## Overview

Product Name	LSM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	LSM1 CASM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>LSM1 homolog, mRNA degradation associated(LSM1) Homo sapiens This gene encodes a member of the LSm family of RNA-binding proteins. LSm proteins form stable heteromers that bind specifically to the 3'-terminal oligo(U) tract of U6 snRNA and may play a role in pre-mRNA splicing by mediating U4/U6 snRNP formation. Increased expression of this gene may play a role in cellular transformation and the progression of several malignancies including lung cancer, mesothelioma and breast cancer. Alternatively spliced transcript variants have been observed for this gene, and a pseudogene of this gene is located on the short arm of chromosome 9. [provided by RefSeq, Nov 2011],</p>