

CHD1 Polyclonal Antibody

Cat No: HR1AP10132

For research use only

Overview

Product Name	CHD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	188kDa
GeneID?Human?	CHD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chromodomain helicase DNA binding protein 1(CHD1) Homo sapiens The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template. [provided by RefSeq, Jul 2008],</p>