

TOPB1 Polyclonal Antibody

Cat No: HR1AP10125

For research use only

Overview

Product Name	TOPB1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	167kDa
GeneID?Human?	TOPBP1 KIAA0259
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>topoisomerase (DNA) II binding protein 1(TOPBP1) Homo sapiens This gene encodes a binding protein which interacts with the C-terminal region of topoisomerase II beta. This interaction suggests a supportive role for this protein in the catalytic reactions of topoisomerase II beta through transient breakages of DNA strands. [provided by RefSeq, Jul 2008],</p>