

RRAS Polyclonal Antibody

Cat No: HR1AP10122

For research use only

Overview

Product Name	RRAS Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RRAS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	related RAS viral (r-ras) oncogene homolog(RRAS) Homo sapiens The protein encoded by this gene is a small GTPase involved in diverse processes including angiogenesis, vascular homeostasis and regeneration, cell adhesion, and neuronal axon guidance. Mutations in this gene are found in many invasive cancers. [provided by RefSeq, Jul 2015],