

## SL9A3 Polyclonal Antibody

Cat No: HR1AP10117

For research use only

### Overview

Product Name	SL9A3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	91kDa
GeneID?Human?	SLC9A3 NHE3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 9 member A3(SLC9A3) Homo sapiens The protein encoded by this gene is an epithelial brush border Na/H exchanger that uses an inward sodium ion gradient to expel acids from the cell. Defects in this gene are a cause of congenital secretory sodium diarrhea. Pseudogenes of this gene exist on chromosomes 10 and 22. [provided by RefSeq, Mar 2016],